

Apc Rat Anti-Mouse Cd45R/B220 Monoclonal Antibody For Flow Cytometry - P06800

Item Number: CM0005861

Introduction

APC-conjugated rat monoclonal antibody targeting mouse CD45R/B220 (leukocyte common antigen, P06800). Validated for flow cytometry (FC) to detect B cells and study T-cell activation. Cross-reacts with mouse.

[Learn More](#)

Parameter	Value
Product Name	APC Rat anti-Mouse CD45R/B220 mAb
Remarks/Alias	Ioc; B220; Cd45; L-CA; Ly-5; T200; CD45R; Lyt-4
Species	Mouse
GeneID (Human)	5788
GeneID	19264
Immunogen	Recombinant protein Recombinant protein of mouse CD45R/B220.
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06800
Protein Weight	130kDa/140kDa/144kDa
Shipping	Ice bag