

Apc Rat Anti-Mouse Cd86 Monoclonal Antibody For Flow Cytometry - P42082

Item Number: CM0005298

Introduction

APC-conjugated rat monoclonal antibody targeting mouse CD86 (B7-2) for flow cytometry. CD86 is a costimulatory molecule on antigen-presenting cells critical for T-cell activation and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	APC Rat anti-Mouse CD86 mAb
Remarks/Alias	B7, B70, MB7, B7-2, B7.2, CLS1, Ly58, ETC-1, Ly-58, MB7-2, Cd28l2, TS/A-2, CD86
Species	Mouse
GeneID (Human)	942
GeneID (Mouse)	12524
Immunogen	Recombinant protein of mouse CD86
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P42082
Protein Weight	35kDa
Shipping	Ice bag