

Apc Rat Anti-Mouse F4/80 Monoclonal Antibody For Flow Cytometry - Q61549

Item Number: CM0005846

Introduction

High-quality APC-conjugated rat monoclonal antibody specifically targeting mouse F4/80 (Adgre1/EMR1), validated for flow cytometry. Ideal for macrophage detection in mouse samples.

[Learn More](#)

Parameter	Value
Product Name	APC Rat anti-Mouse F4/80 mAb
Remarks/Alias	Emr1; Ly71; F4/80; Gp480; TM7LN3; DD7A5-7; EGF-TM7; Adgre1
Species	Mouse
GenelD (Human)	2015
GenelD	13733
Immunogen	Recombinant protein of mouse F4/80.
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q61549
Protein Weight	102kDa
Shipping	Ice bag