

Apc Rat Anti-Mouse Ly-6G/Ly-6C (Gr-1) Monoclonal Antibody For Fc - P35461 / P0Cw02

Item Number: CM0007912

Introduction

APC-conjugated rat monoclonal antibody targeting mouse Ly-6G/Ly-6C (Gr-1) validated for flow cytometry (FC). Suitable for immunophenotyping of granulocytes and monocytes in mouse bone marrow and peripheral blood samples.

[Learn More](#)

Parameter	Value
Product Name	APC Rat anti-Mouse Ly-6G/Ly-6C (Gr-1) mAb
Remarks/Alias	Ly6c; Ly-6C; Ly-6C1
Species	Mouse
GenelD	546644 (Ly6g), 17067 (Ly6c)
Immunogen	A synthesized peptide derived from mouse Ly-6G/Ly-6C (Gr-1)
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P35461 / P0CW02
Protein Weight	15kDa / 8kDa
Shipping	Ice bag