

Apc Rat Anti-Mouse Ly6G Monoclonal Antibody For Flow Cytometry - P35461

Item Number: CM0008604

Introduction

APC-conjugated rat anti-mouse Ly6G monoclonal antibody, validated for flow cytometry (FC). Reacts with mouse Ly6G, a GPI-anchored cell surface marker. Suitable for immunophenotyping and bone marrow studies.

[Learn More](#)

Parameter	Value
Product Name	APC Rat anti-Mouse Ly6G mAb
Remarks/Alias	Gr1; Gr-1; Ly-6G
Species	Mouse
GenelD (human)	80740
GenelD	546644
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-105 of mouse Ly-6G (NP_001297367.1)
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P35461
Protein Weight	14kDa
Shipping	Ice bag