

Apc Rat Anti-Mouse Ter-119 Monoclonal Antibody For Flow Cytometry - Ter-119

Item Number: CM0012626



Introduction

This APC Rat anti-Mouse TER-119 monoclonal antibody is designed for flow cytometry detection of mouse TER-119 antigen on erythroid cells. Ideal for immunophenotyping and erythroid development studies.

[Learn More](#)

Parameter	Value
Product Name	APC Rat anti-Mouse TER-119 mAb
Remarks/Alias	Ter119, TER-119
Species	Mouse
Gene ID	104231
Immunogen	Recombinant protein of mouse TER-119
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
Shipping	Ice bag