

Apc Streptavidin Monoclonal Antibody For Flow Cytometry - Apc Streptavidin

Item Number: CM0004592

Introduction

APC Streptavidin (SA V1, SA V2) monoclonal antibody conjugate designed for flow cytometry applications. Species-independent cross reactivity ensures broad compatibility. Reliable raw material for IVD assay development.

[Learn More](#)

Parameter	Specification
Product Name	APC Streptavidin
Remarks/Alias	SA V1, SA V2
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Species independent
Shipping	Ice bag