

Apc Streptavidin

Item Number: CM0011197

Introduction

APC Streptavidin is a rabbit-derived streptavidin conjugate designed for flow cytometry (FC) applications with species-independent cross reactivity. Ideal for detecting biotinylated probes, it supports reliable IVD and research workflows.

[Learn More](#)

Parameter	Value
Product Name	APC Streptavidin
Remarks/Alias	SA V1, SA V2
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Species independent
Shipping	Ice bag