


Apc/Cyanine7 Anti-Axl Monoclonal Antibody For Flow Cytometry - P30530

Item Number: CM0008169



Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody targeting human AXL receptor tyrosine kinase. Validated for flow cytometry, suitable for detecting AXL in oncology and virology research. 98 kDa, cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Human AXL mAb
Remarks/Alias	ARK; UFO; JTK11; Tyro7
Species	Human
GeneID (Human)	558
GeneID	558
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 33-449 of human AXL (NP_068713.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P30530
Protein Weight	98kDa
Shipping	Ice bag