

Apc/Cyanine7 Anti-Cd185/Cxcr5 Monoclonal Antibody For Fc - P32302

Item Number: CM0005676

Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody targeting human CD185/CXCR5 (Burkitt lymphoma receptor 1). Validated for flow cytometry, this antibody enables detection of CXCR5, a key B-cell migration receptor, in human samples.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Human CD185/CXCR5 mAb
Remarks/Alias	BLR1; CD185; MDR15
Species	Human
GeneID (Human)	643
GeneID	643
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-47 of human CD185/CXCR5 (NP_033981.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P32302
Protein Weight	42 kDa
Shipping	Ice bag