

Apc/Cyanine7 Anti-Cd20 Monoclonal Antibody For Flow Cytometry - P11836

Item Number: CM0012456



Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody targeting human and cynomolgus CD20, validated for flow cytometry. Ideal for B-cell immunophenotyping and diagnostic research.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Human/Monkey CD20 mAb
Remarks/Alias	B1, S7, Bp35, CD20, FMC7, CVID5, LEU-16
Species	Human
Gene ID (Human)	931
Gene ID	931
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human CD20 mAb (NP_068769.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P11836
Protein Weight	14kDa/33kDa
Shipping	Ice bag