

APC/Cyanine7 Anti-Cd3E Rabbit Monoclonal Antibody For Intracellular Flow Cytometry - P07766

Item Number: CM0004230



Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody targeting human CD3 epsilon chain. Validated for intracellular flow cytometry (FC intra); ideal for T-cell immunophenotyping and TCR signaling research.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Human cCD3 mAb
Remarks/Alias	T3E; TCRE; IMD18; CD3epsilon
Species	Human
GeneID (Human)	916
GeneID	916
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human CD3E (NP_000724.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P07766
Protein Weight	23kDa
Shipping	Ice bag