

Apc/Cyanine7 Anti-Cd45 Monoclonal Antibody For Flow Cytometry - P08575-8

Item Number: CM0011992

Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody targeting isoform 6 of human/monkey CD45 (PTPRC). Validated for flow cytometry and cross-reactive with human and cynomolgus species.

[Learn More](#)

Product Name	APC/Cyanine7 Rabbit anti-Human/Monkey CD45 mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
GeneID (Human)	5788
GeneID	5788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-577 of human CD45 (P08575)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P08575-8
Protein Weight	138kDa
Shipping	Ice bag