

APC/Cyanine7 Anti-Hla-Dra Monoclonal Antibody For Flow Cytometry - P01903

Item Number: CM0008845



Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody targeting human HLA-DRA, validated for flow cytometry on human and cynomolgus monkey cells. Ideal for MHC class II immunophenotyping and T cell response studies.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Human/Monkey HLA-DR mAb
Remarks/Alias	HLA-DRA1
Species	Human, Cynomolgus monkey
GeneID (Human)	3122
GeneID	3122
Immunogen	Recombinant protein of human HLA-DR
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01903
Protein Weight	29kDa
Shipping	Ice bag