

APC/Cyanine7 Anti-Human Cd38 Rabbit Monoclonal Antibody For Flow Cytometry - P28907

Item Number: CM0010588



Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody targeting human CD38 (ADPRC1/cADPR1). Suitable for flow cytometry (FC) applications, validated for human samples. Ideal for immunology, cancer, and diabetes research.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Human CD38 mAb
Remarks/Alias	ADPRC1, cADPR1, ADPRC 1
Species	Human
GeneID (Human)	952
GeneID	952
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 43-300 of human CD38 (P28907)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P28907
Protein Weight	13kDa/34kDa
Shipping	Ice bag