

APC/Cyanine7 Anti-Human/Monkey Cd127/IL-7R α Monoclonal Antibody For Flow Cytometry - P16871

Item Number: CM0004319



Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody targeting human and cynomolgus CD127/IL-7 receptor alpha. Validated for flow cytometry, ideal for immunology and lymphocyte development studies.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Human/Monkey CD127/IL-7R α mAb
Remarks/Alias	ILRA; CD127; IL7RA; CDW127; IMD104; sIL-7R; Inc-IL7R; IL7Ralpha; IL-7Ralpha; IL-7R-alpha
Species	Human
Gene ID	3575
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-236 human CD127/IL-7R α (NP_002176.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P16871
Protein Weight	52kDa
Shipping	Ice bag