

Apc/Cyanine7 Anti-Mouse Cd19 Monoclonal Antibody For Flow Cytometry - P25918

Item Number: CM0011195

Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody targeting mouse CD19 for flow cytometry. Detects B-lymphocyte antigen on B cells; supports immunology and B-cell research. - P25918.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Mouse CD19 mAb
Remarks/Alias	Cd19
Species	Mouse
GenelD (Human)	930
GenelD	12478
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 19-287 of mouse CD19 (NP_033974.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P25918
Protein Weight	60kDa
Shipping	Ice bag