

APC/Cyanine7 Anti-Mouse Cd38 Monoclonal Antibody For Flow Cytometry - P56528

Item Number: CM0009597



Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody against mouse CD38, validated for flow cytometry. Suitable for immunology research involving lymphocyte activation, calcium signaling, and cell adhesion. - P56528

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Mouse CD38 mAb
Remarks/Alias	I-19; ADPRC 1; Cd38-rs1
Species	Mouse
GeneID (Human)	952
GeneID	12494
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P56528
Protein Weight	34kDa
Shipping	Ice bag