

Apc/Cyanine7 Anti-Mouse Cd45 Rabbit Monoclonal Antibody For Flow Cytometry - P06800

Item Number: CM0008390

Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody targeting mouse CD45 (PTPRC, leukocyte common antigen, Ly-5, B220). Validated for flow cytometry applications in mouse samples. Ideal for immunophenotyping and immune cell analysis.

[Learn More](#)

Parameter	Specification
Product Name	APC/Cyanine7 Rabbit anti-Mouse CD45 mAb
Remarks/Alias	loc; B220; Cd45; L-CA; Ly-5; T200; CD45R; Lyt-4
Species	Mouse
Gene ID (Human)	5788
Gene ID	19264
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-566 of mouse CD45 (NP_001104786.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06800
Protein Weight	130kDa/140kDa/144kDa
Shipping	Ice bag