

Apc/Cyanine7 Anti-Mouse Cd69 Monoclonal Antibody For Flow Cytometry - P37217

Item Number: CM0007941

Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody targeting mouse CD69 (Early activation antigen CD69). Validated for flow cytometry (FC). Suitable for immune cell activation studies and T-cell research. Reacts with mouse. Molecular weight ~22kDa.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Mouse CD69 mAb
Remarks/Alias	AIM; VEA; 5830438K24Rik
Species	Mouse
GeneID (human)	969
GeneID	12515
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to aminoacids 64-199 of mouse CD69(NP_001028294.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P37217
Protein Weight	22kDa
Shipping	Ice bag