

APC/Cyanine7 Mouse Anti-Human Cd3 Monoclonal Antibody For Fc - P07766

Item Number: CM0007870



Introduction

APC/Cyanine7-conjugated mouse monoclonal antibody targeting human CD3 epsilon chain for flow cytometry. Reacts with human CD3, part of the TCR-CD3 complex. Ideal for T-cell subset analysis and immunological research.

[Learn More](#)

Specification	Details
Product Name	APC/Cyanine7 Mouse anti-Human CD3 mAb
Remarks/Alias	T3E; TCRE; IMD18; CD3epsilon
Species	Human
GenelD	916
Immunogen	Recombinant protein of human CD3
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P07766
Protein Weight	23kDa
Shipping	Ice bag