

APC/Cyanine7 Rabbit Anti-Human Cd10 Monoclonal Antibody For Fc - P08473

Item Number: CM0005346



Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody against human CD10 (Neprilysin, CALLA), validated for flow cytometry. Ideal for studying enkephalin degradation, cardiovascular signaling, and neurodegenerative diseases.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Human CD10 mAb
Remarks/Alias	NEP; SFE; CD10; CALLA; CMT2T; SCA43
Species	Human
GeneID (Human)	4311
GeneID	4311
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 52-750 of human CD10 (NP_000893.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08473
Protein Weight	86 kDa
Shipping	Ice bag