

APC/Cyanine7 Rabbit Anti-Human Cd56/Ncam1 Monoclonal Antibody For Flow Cytometry - P13591

Item Number: CM0005874



Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody targeting human CD56/NCAM1 specifically for flow cytometry applications. Recognizes the neural cell adhesion molecule 1 (NCAM1), a ~95 kDa membrane glycoprotein mediating cell adhesion and serving as a viral entry receptor. Designed for human samples.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Human CD56/NCAM1 mAb
Remarks/Alias	CD56; NCAM; MSK39
Species	Human
GeneID (Human)	4684
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-413 of human CD56/NCAM1 (NP_851996.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P13591
Protein Weight	95kDa
Shipping	Ice bag