

APC/Cyanine7 Rabbit Anti-Human CD61 Mab For Flow Cytometry - P05106

Item Number: CM0008544



Introduction

APC/Cyanine7-conjugated rabbit anti-human CD61 (integrin beta-3, GPIIIa) monoclonal antibody for flow cytometry. 87 kDa, human cross-reactivity. Supports platelet aggregation, integrin signaling, and cardiovascular disease research.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Human CD61 mAb
Remarks/Alias	GT; GT2; CD61; GP3A; BDPLT2; GPIIIa; BDPLT16; BDPLT24
Species	Human
Gene ID (Human)	3690
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-718 of human CD61 (NP_000203.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P05106
Protein Weight	87 kDa
Shipping	Ice bag