


Apc/Cyanine7 Rabbit Anti-Human Cd71 Monoclonal Antibody For Fc - P02786

Item Number: CM0011293



Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody targeting human CD71 (TFRC, T9, TR, TFR, p90, CD71, TFR1, TRFR). Suitable for flow cytometry detection of transferrin receptor 1 in human samples. Molecular weight approximately 85kDa.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Human CD71 mAb
Remarks/Alias	T9; TR; TFR; p90; CD71; TFR1; TRFR; IMD46
Species	Human
GeneID (Human)	7037
GeneID	7037
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-250 of human CD71 (NP_003225.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P02786
Protein Weight	85kDa
Shipping	Ice bag