


Apc/Cyanine7 Rabbit Anti-Human Egfr Monoclonal Antibody For Fc - P00533

Item Number: CM0009425



Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody against human EGFR (HER1/ErbB-1). Validated for flow cytometry, enabling precise detection of cell surface EGFR in oncology and signaling studies.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Human EGFR mAb
Remarks/Alias	ERBB; ERRP; HER1; mENA; ERBB1; PIG61; NISBD2
Species	Human
GeneID (Human)	1956
GeneID	1956
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 335-525 of human EGFR(NP_005219.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P00533
Protein Weight	44kDa/69kDa/77kDa/134kDa
Shipping	Ice bag