

APC/Cyanine7 Rabbit Anti-Human Ifn-Gamma Monoclonal Antibody For Flow Cytometry (Intracellular) - P01579

Item Number: CM0008129



Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody targeting human IFN-gamma. Suitable for intracellular flow cytometry, it aids research in immune response, antiviral defense, and immunotherapy. Validated with immunogen spanning amino acids 24-166.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Human IFN-gamma mAb
Remarks/Alias	IFG; IFI; IMD69
Species	Human
Gene ID (Human)	3458
Gene ID	3458
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 24-166 of human IFN-gamma (NP_000610.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P01579
Protein Weight	19kDa
Shipping	Ice bag