

APC/Cyanine7 Rabbit Anti-Human/Monkey Cd14 Monoclonal Antibody For Fc - P08571

Item Number: CM0009023



Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody targeting CD14 (human/monkey). Validated for flow cytometry (FC). Ideal for studying innate immunity, LPS signaling, and monocyte/macrophage biology.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Human/Monkey CD14 mAb
Remarks/Alias	CD14
Species	Human
GeneID (human)	929
GeneID	929
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P08571
Protein Weight	40kDa
Shipping	Ice bag