

APC/Cyanine7 Rabbit Anti-Human/Monkey Cd19 Mab For Flow Cytometry - Cd19

Item Number: CM0008488

Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody targeting human and cynomolgus monkey CD19. Specifically validated for flow cytometry (FC) applications. Ideal for B-cell immunophenotyping in immunology and diagnostic research. Cross-reactive with human and cynomolgus monkey.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Human/Monkey CD19 mAb
Species	Human, Monkey
Gene ID (Human)	930
Gene ID (Cynomolgus Monkey)	102145514
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-292 of monkey CD19.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
Protein Weight	63kDa
Shipping	Ice bag