

APC/Cyanine7 Rabbit Anti-Human/Monkey Cd25 Monoclonal Antibody For Fc - P01589

Item Number: CM0009288



Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody targeting human and cynomolgus CD25 (IL-2 receptor alpha, p55). Designed for flow cytometry applications in immune regulation studies, including regulatory T cell (Treg) research. Lot-specific validation recommended.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Human/Monkey CD25 mAb
Remarks/Alias	p55; CD25; IL2R; IMD41; TCGFR; IDDM10
Species	Human
GeneID (human)	3559
GeneID	3559
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-213 of human CD25 (NP_000408.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01589
Protein Weight	31kDa
Shipping	Ice bag