


APC/Cyanine7 Rabbit Anti-Human/Monkey Cd38 Polymab® - P28907

Item Number: CM0008964



Introduction

APC/Cyanine7-conjugated rabbit anti-human/monkey CD38 monoclonal antibody (PolymAb®) for flow cytometry. Detects CD38 in human and cynomolgus models. Ideal for immune cell profiling and CD38 research.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Human/Monkey CD38 PolymAb®
Remarks/Alias	ADPRC1; cADPR1; ADPRC 1
Species	Human
GeneID (Human)	952
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 43-300 of human CD38 (P28907)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P28907
Protein Weight	34kDa
Shipping	Ice bag