


APC/Cyanine7 Rabbit Anti-Human/Monkey Cd4 Monoclonal Antibody For Fc - P01730

Item Number: CM0012198



Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody targeting human and cynomolgus CD4, validated for flow cytometry. Enables identification of T-helper cells and immune research applications.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Human/Monkey CD4 mAb
Remarks/Alias	T4; IMD79; Leu-3; OKT4D; CD4mut
Species	Human
GeneID (Human)	920
GeneID	920
Immunogen	Recombinant protein of human CD4
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01730
Protein Weight	51kDa
Shipping	Ice bag