


APC/Cyanine7 Rabbit Anti-Human/Monkey Cd8A Monoclonal Antibody For Flow Cytometry - P01732

Item Number: CM0012747



Introduction

APC/Cyanine7-conjugated rabbit anti-CD8a monoclonal antibody validated for flow cytometry. Recognizes human and cynomolgus CD8a; supports T-cell and immune response research.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Human/Monkey CD8a mAb
Remarks/Alias	CD8; p32; Leu2; CD8alpha
Target Species	Human
Gene ID (Human)	925
Gene ID	925
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-182 of human CD8a (NP_001759.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01732
Protein Weight	21kDa/25kDa/30kDa
Shipping	Ice bag