


APC/Cyanine7 Rabbit Anti-Human/Monkey IgD Monoclonal Antibody For Flow Cytometry - P01880

Item Number: CM0008027



Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody targeting human and cynomolgus monkey immunoglobulin delta chain for flow cytometry. Ideal for B-cell immunophenotyping and humoral immunity studies.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Human/Monkey IgD mAb
Remarks/Alias	Immunoglobulin heavy constant delta
Species	Human
GeneID (Human)	3495
GeneID	3495
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 99-384 of human IgD.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01880
Protein Weight	47kDa
Shipping	Ice bag