


APC/Cyanine7 Rabbit Anti-Mouse Cd11B Monoclonal Antibody For Flow Cytometry - P05555

Item Number: CM0010218



Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody against mouse CD11b (integrin alpha-M) for flow cytometry. Recognizes mouse CD11b, a leukocyte adhesion receptor. Ideal for immunophenotyping and immune cell research.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Mouse CD11b mAb
Remarks/Alias	CR3; CR3A; MAC1; Cd11b; Ly-40; Mac-1; Mac-1a; CD11b/CD18; F730045J24Rik
Species	Mouse
GeneID (Human)	3684
GeneID	16409
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 143-337 of mouse CD11b (NP_032427.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P05555
Protein Weight	127kDa
Shipping	Ice bag