

APC/Cyanine7 Rabbit Anti-Mouse CD19 Mab For Fc - P25918

Item Number: CM0008308

Introduction

This APC/Cyanine7-conjugated rabbit anti-mouse CD19 monoclonal antibody is intended for flow cytometric detection of mouse CD19, a B-cell lineage marker. Ideal for immunophenotyping and immunological studies.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Mouse CD19 mAb
Species	Mouse
Gene ID (Human)	930
Gene ID (Mouse)	12478
Immunogen	Recombinant protein of mouse CD19.
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry (FC)
Cross Reactivity	Mouse
SWISS	P25918
Protein Weight	60kDa
Shipping	Ice bag