

Apc/Cyanine7 Rabbit Anti-Mouse Cd25 Monoclonal Antibody For Flow Cytometry - P01590

Item Number: CM0004646

Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody targeting mouse CD25 (IL-2RA, p55). Validated for flow cytometry. This reagent enables precise detection of cell surface CD25 on murine immune cells, supporting research in immune tolerance and regulatory T cell (Treg) biology.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Mouse CD25 mAb
Remarks/Alias	CD25; IL2r; Ly-43
Species	Mouse
GeneID (human)	3559
GeneID	16184
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-236 of mouse CD25 (NP_032393.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01590
Protein Weight	31kDa
Shipping	Ice bag