

APC/Cyanine7 Rabbit Anti-Mouse Cd44 Monoclonal Antibody For Flow Cytometry - P15379

Item Number: CM0004622



Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody against mouse CD44 (Ly-24, Pgp-1, HERMES). Validated for flow cytometry. Reacts with mouse. Predicted molecular weight 85 kDa. UniProt P15379.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Mouse CD44 mAb
Remarks/Alias	Ly-24; Pgp-1; HERMES
Species	Mouse
GeneID (Human)	960
GeneID	12505
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-176 of mouse CD44 (NP_033981.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P15379
Protein Weight	85 kDa
Shipping	Ice bag