

APC/Cyanine7 Rabbit Anti-Mouse CD8A Monoclonal Antibody For Flow Cytometry - P01731

Item Number: CM0005535



Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody against mouse CD8a (Ly-2/Lyt-2) for flow cytometry. Recognizes the 27 kDa T-cell coreceptor; validated for mouse CD8a detection in FC.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-Mouse CD8a mAb
Remarks/Alias	Ly-2; Ly-B; Ly-35; Lyt-2
Species	Mouse
GeneID (Human)	941
GeneID	12525
Immunogen	Recombinant protein containing a sequence corresponding to amino acids 28-196 of mouse CD8a (NP_001074579.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01731
Protein Weight	27kDa
Shipping	Ice bag