

Apc/Cyanine7 Rabbit Anti-Mouse Cd90.2 Monoclonal Antibody For Flow Cytometry - P01831

Item Number: CM0007782

Introduction

APC/Cyanine7-conjugated rabbit monoclonal antibody targeting mouse CD90.2 (Thy-1), validated for flow cytometry. GPI-anchored membrane glycoprotein involved in cell interactions. Ideal for immunological and neuroscience research.

[Learn More](#)

Parameter	Specification
Product Name	APC/Cyanine7 Rabbit anti-Mouse CD90.2 mAb
Remarks/Alias	T25; CD90; Thy-1; Thy1.1; Thy1.2; Thy-1.2
Species	Mouse
Gene ID (Human)	7070
Gene ID	21838
Immunogen	Synthetic peptide: A synthesized peptide derived from mouse CD90.2
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01831
Protein Weight	18kDa
Shipping	Ice bag