

Apc/Cyanine7 Rabbit Anti-Nhp Igd Monoclonal Antibody - Cynomolgus Monkey Igd

Item Number: CM0005941

Introduction

This APC/Cyanine7-conjugated rabbit monoclonal antibody targets non-human primate (NHP) immunoglobulin delta heavy chain and is intended for flow cytometry (FC). It cross-reacts with cynomolgus and rhesus monkey IgD, supporting immune cell analysis in primate research models.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rabbit anti-NHP IgD mAb
Species	Cynomolgus monkey
Immunogen	Recombinant protein of cynomolgus monkey immunoglobulin delta heavy chain
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Cynomolgus, Rhesus
Shipping	Ice bag