


APC/Cyanine7 Rat Anti-Mouse Mhc Class Ii (I-A/I-E) Monoclonal Antibody For Fc - P18468 / P06342 / Q3U060 / P14435

Item Number: CM0010101



Introduction

APC/Cyanine7-conjugated rat monoclonal antibody targeting mouse MHC Class II (I-A/I-E) for flow cytometry applications in immunological research and diagnostic assay development. Reacts with multiple haplotypes.

[Learn More](#)

Parameter	Value
Product Name	APC/Cyanine7 Rat anti-Mouse MHC Class II (I-A/I-E) mAb
Remarks/Alias	Eb; Ia4; H2Eb; Ia-4; H-2Eb; IAb; Ia2; Ia-2; Abeta; H-2Ab; H2-Ab; Rmcs1; I-Abeta; Ia5; Ia-5; H-2Eb2; Iebeta2; A130038H09Rik; Ia1; H2Aa; Ia-1; H-2Aa; Aalpha; IAalpha; I-Aalpha
Species	Mouse
GeneID	14969, 14961, 381091, 14960
Immunogen	Recombinant protein of mouse MHC Class II (I-A/I-E)
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P18468 / P06342 / Q3U060 / P14435
Protein Weight	30kDa / 32kDa / 26kDa
Shipping	Ice bag