


Apex1/Ape1 Rabbit Mab For Wb, Ihc-P, Elisa - P27695

Item Number: CM0002932



Introduction

Rabbit monoclonal anti-APEX1/APE1 antibody, validated for WB, IHC-P, and ELISA, with reactivity to human, mouse, and rat. Suitable for DNA repair and oxidative stress studies.

[Learn More](#)

Parameter	Value
Product Name	APEX1/APE1 Rabbit mAb
Remarks/Alias	APE; APX; APE1; APEN; APEX; HAP1; REF1; APEX1/APE1
Species	Human
GeneID (Human)	328
GeneID	328
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human APEX1/APE1 (P27695).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27695
Protein Weight	36kDa
Shipping	Ice bag