


Apoa1 Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - P02647

Item Number: CM0013002



Introduction

Rabbit polyclonal antibody targeting human APOA1 (Apolipoprotein A-I). Validated for WB, IHC-P, IF/ICC, ELISA with cross-reactivity to mouse and rat. Ideal for reverse cholesterol transport and cardiovascular disease studies.

[Learn More](#)

Parameter	Value
Product Name	APOA1 Rabbit pAb
Remarks/Alias	HPALP2; apo(a); APOA1
Species	Human
GeneID (Human)	335
GeneID	335
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-267 of human APOA1 (NP_000030.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02647
Protein Weight	28kDa
Shipping	Ice bag