

Apolipoprotein Ci Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P02654

Item Number: CM0008449



Introduction

Rabbit monoclonal antibody against human Apolipoprotein Ci (APOC1), validated for WB, IF-P, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for lipid metabolism and cardiovascular disease research.

[Learn More](#)

Parameter	Value
Product Name	Apolipoprotein Ci Rabbit mAb
Remarks/Alias	Apo-CI; ApoC-I; apo-CIB; apoC-IB; Apolipoprotein Ci
Species	Human
Gene ID (Human)	341
Gene ID	341
Immunogen	Synthetic Peptide: A synthesized peptide derived from human APOC1(P02654)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02654
Protein Weight	9kDa
Shipping	Ice bag