

Aprr3 Rabbit Polyclonal Antibody For Arabidopsis Thaliana (Wb, Elisa) - Aprr3 Rabbit Pab

Item Number: CM0019236

Introduction

APRR3 rabbit polyclonal antibody raised against recombinant APRR3 protein of Arabidopsis thaliana. Validated for Western blotting and ELISA applications in plant science research.

[Learn More](#)

Parameter	Value
Product Name	APRR3 Rabbit pAb
Remarks/Alias	APRR3; MGO3.8; MGO3_8; pseudo-response regulator 3
Species	Arabidopsis thaliana
Gene ID	836132
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 245-371 of Arabidopsis thaliana APRR3 (NP_568919.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Arabidopsis thaliana
Shipping	Ice bag