

Aprr7 Rabbit Polyclonal Antibody - Aprr7

Item Number: CM0019936

Introduction

High-quality APRR7 Rabbit Polyclonal Antibody validated for WB and ELISA. Specifically detects Arabidopsis thaliana APRR7 (pseudo-response regulator 7). Ideal for plant signaling research.

[Learn More](#)

Parameter	Value
Product Name	APRR7 Rabbit pAb
Remarks/Alias	APRR7; F9G14.120; F9G14_120; pseudo-response regulator 7
Species	Arabidopsis thaliana
Gene ID	831793
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 236-339 of Arabidopsis thaliana APRR7 (NP_568107.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Arabidopsis thaliana
Shipping	Ice bag