


Aqp5 Rabbit Pab For Wb, Ihc-P, Elisa - P55064

Item Number: CM0014877



Introduction

CamelBio offers AQP5 Rabbit pAb, a polyclonal antibody targeting human Aquaporin-5 (P55064) for WB, IHC-P, ELISA, with cross-reactivity to mouse and rat. Ideal for water transport studies.

[Learn More](#)

Parameter	Value
Product Name	AQP5 Rabbit pAb
Remarks/Alias	PPKB; AQP-5; AQP5
Species	Human
Gene ID (Human)	362
Gene ID	362
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-265 of human AQP5 (NP_001642.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P55064
Protein Weight	28kDa
Shipping	Ice bag