

Ar-V7 Rabbit Pab For Wb, If/Icc, Elisa - P10275-3

Item Number: CM0023357

Introduction

Rabbit polyclonal antibody against human AR-V7 (androgen receptor splice variant) for WB, IF/ICC, ELISA. Targets constitutively active isoform lacking ligand-binding domain, 67 kDa. Validated in human samples.

[Learn More](#)

Parameter	Value
Product Name	AR-V7 Rabbit pAb
Remarks/Alias	KD; AIS; AR8; TFM; DHTR; SBMA; HYPSP1; NR3C4; SMAX1; HUMARA; AR-V7
Species	Human
GeneID (Human)	367
GeneID	367
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 544-644 of human AR-V7 (NP_001334990.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P10275-3
Protein Weight	67kDa
Shipping	Ice bag