

Arginase-1 Rabbit Mab For Wb, If-P, Ihc-P, Elisa - Q61176

Item Number: CM0007929

Introduction

Rabbit monoclonal antibody targeting mouse Arginase-1 (Arg-1), validated for WB, IF-P, IHC-P, and ELISA. Reacts with human, mouse, and rat. Ideal for studying urea cycle, immune regulation, and tumor immunity.

[Learn More](#)

Parameter	Value
Product Name	Arginase-1 Rabbit mAb
Remarks/Alias	AI; PGIF; Arg-1; Arginase-1
Species	Mouse
GeneID (Human)	383
GeneID	11846
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of mouse Arginase-1 (NP_031508.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q61176
Protein Weight	35kDa
Shipping	Ice bag